




KONICA MINOLTA

Traceability System for Spectroradiometer

January, 2004

Konica Minolta Photo Imaging Europe GmbH
Instrument Systems Division



Interval between calibration:

Luminous Intensity / Illuminance Standard Lamp: Each time a product is calibrated

Diffuse Transmission Plate: Each time a product is calibrated

Standard White Plate: 5 years

Wavelength Calibration Filter: 8 years

Spectroscope: 1 year

JCSS0026: Japan Calibration Service System, the calibration laboratory accreditation system in Japan.

The organization number for KONICA MINOLTA SENSING, INC. is 0026